



Email : sales@laboratoryinstrumentindia.com

Product Name :
A Class Petrol Engine With Injection

Product Code :
LBNY-0005-12000078



Description :

The start can be engaged, the transmission shifted. The engine is driven by a 220V geared motor, all the assemblies also run.

The following cutaways are:
cylinder head, cylinder block, crankcase, oil pan, one piston, roller valve lever with clearance compensation element, generator (internally ventilated), starter, oil pump, chain case, intake manifold, electronic control unit, actuator, air filter box, double ignition coil, throttle valve, thermostat, air-mass sensor, valve cover, oil filter, 5-gear transmission and differential.

Technical Specification :

A Class Petrol Engine With Injection



Laboratory Instrument India