




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
High Pressure Pump For Engines With Direct Petrol Injection

**Product Code :**  
LBNY-0005-12000086



**Description :**

The eccentric shaft is guided on the drive side by a roller bearing and opposite by a plastic friction bearing. Also cutaway: fuel return with shut-off attenuator, quantity control valve with magnetic coil and filter. The necessary high pressure is generated by a 3-plunger radial-piston pump with the following features: compact construction, stainless steel housing with integrated quantity control valve, fuel-lubricated, shut-off attenuator, max. fuel pressure of 200 bar. The knurled screw can be used to drive the high-pressure pump and one sees how the pump pistons work.

**Technical Specification :**

High Pressure Pump For Engines With Direct Petrol Injection



**Laboratory Instrument India**