




Email : sales@laboratoryinstrumentindia.com

Product Name :
Exhaust Gas Turbocharger

Product Code :
LBNY-0005-12000090



Description :

The two sides are marked by colours: compression side - blue, turbine side - red. All parts can be easily moved. The housing of the charger is cut in a way you can see compression wheel, turbine wheel and bearings. The housing of the charger is cut in a way you can see compression wheel, turbine wheel and bearings.

Technical Specification :

Exhaust Gas Turbocharger



Laboratory Instrument India