




Email : sales@laboratoryinstrumentindia.com

Product Name :
Diesel Powered Industrial Motor

Product Code :
LBNY-0005-1200097



Description :

The engine can be turned easily. This engine is a 1-cylinder naturally aspirated diesel engine with balancer shaft and bottom-mounted camshaft.

Cutaway are: engine case, injection nozzle, injection pump, cylinder, cylinder head, muffler, starter, air filter and rotor.

Technical Specification :

Diesel Powered Industrial Motor



Laboratory Instrument India