




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Demand Driven Oil Pump With Suction Stage

**Product Code :**  
LBNY-0005-120000100



**Description :**

This control concept allows the engine lubricating and cooling points to be supplied independently of engine load and engine speed absolutely and with a significantly lower drive power particularly in the part-load range than would be possible with a non. A newly developed vane-type oil pump with demand-driven quantity regulation and two map-controlled pressure stages is used in the new 6- and 8-cylinder V engines. Controlled pump. The pump is easy to turn.

The following are easy to recognise:

- Valve for motor oil pump
- Rotor with vane
- Control piston
- Cold-start protection valve
- Toothed wheel for oil pump
- Adjusting ring.

**Technical Specification :**

Demand Driven Oil Pump With Suction Stage

---

# Laboratory Instrument India

