



Email : sales@laboratoryinstrumentindia.com

Product Name :
Camshaft Adjustment

Product Code :
LBNY-0005-120000102



Description :

The hydraulic piston and the piston of the solenoid valve are cut away
Drive of the inlet camshaft via a control chain
Drive of the outlet camshaft via a rotary button
Holding the outlet camshaft with the help of a set screw
The function of the chain tensor is easily recognised
With the help of a push bar, the hydraulic piston can be displaced and the adjustment of the inlet camshaft read off directly on a scale.

Technical Specification :

Camshaft Adjustment



Laboratory Instrument India