



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Camshaft Adjustment

**Product Code :**  
LBNY-0005-120000102



**Description :**

The hydraulic piston and the piston of the solenoid valve are cut away  
Drive of the inlet camshaft via a control chain  
Drive of the outlet camshaft via a rotary button  
Holding the outlet camshaft with the help of a set screw  
The function of the chain tensor is easily recognised  
With the help of a push bar, the hydraulic piston can be displaced and the adjustment of the inlet camshaft read off directly on a scale.

**Technical Specification :**

Camshaft Adjustment



**Laboratory Instrument India**