



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Common Rail High Pressure Pump With Fuel Delivery  
Regulation Valve

**Product Code :**  
LBNY-0005-120000103



**Description :**

The fuel-delivery regulation valve and the fuel channels are also cutaway  
This high-pressure pump is the latest generation  
By turning the star grip, one sees the function of the pump plungers  
These pumps are installed in the latest Daimler-Benz-CDI vehicles  
Pump cylinder and housing from a block  
Injection pressure: 1500 bar

**Technical Specification :**

Common Rail High Pressure Pump With Fuel Delivery Regulation Valve



**Laboratory Instrument India**