



Email : sales@laboratoryinstrumentindia.com

Product Name :
Common Rail High Pressure Pump With Fuel Delivery
Regulation Valve

Product Code :
LBNY-0005-120000103



Description :

The fuel-delivery regulation valve and the fuel channels are also cutaway
This high-pressure pump is the latest generation
By turning the star grip, one sees the function of the pump plungers
These pumps are installed in the latest Daimler-Benz-CDI vehicles
Pump cylinder and housing from a block
Injection pressure: 1500 bar

Technical Specification :

Common Rail High Pressure Pump With Fuel Delivery Regulation Valve



Laboratory Instrument India