



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Air Flow Sensor Platinum Heating Filament

**Product Code :**  
LBNY-0005-120000104



**Description :**

The printed circuit board with its hybrid circuit can be seen at the back on the cut-away housing. Inside the cut-away plastic housing, the platinum heating filament, the precision resistor and the temperature compensation resistor can be seen in the inner tube.

**Technical Specification :**

Air Flow Sensor Platinum Heating Filament



**Laboratory Instrument India**