



Email : sales@laboratoryinstrumentindia.com

Product Name :
Common Rail High Pressure Pump With Fuel Delivery
Regulation Valve And Injector

Product Code :
LBNY-0005-120000107



Description :

The fuel-delivery regulation valve and the fuel channels are also cutaway.
This high-pressure pump is the latest generation
By turning the star grip, one sees the function of the pump plungers
These pumps are installed in the latest vehicles.
Pump cylinder and housing from a block
The following can easily be seen in the injector: solenoid valve, high and low-pressure channels, injection nozzle, feed and drain bores.

Technical Specification :

Common Rail High Pressure Pump With Fuel Delivery Regulation Valve And Injector



Laboratory Instrument India