



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Common Rail High Pressure Pump With Fuel Delivery  
Regulation Valve And Injector

**Product Code :**  
LBNY-0005-120000107



**Description :**

The fuel-delivery regulation valve and the fuel channels are also cutaway.  
This high-pressure pump is the latest generation  
By turning the star grip, one sees the function of the pump plungers  
These pumps are installed in the latest vehicles.  
Pump cylinder and housing from a block  
The following can easily be seen in the injector: solenoid valve, high and low-pressure channels, injection nozzle, feed and drain bores.

**Technical Specification :**

Common Rail High Pressure Pump With Fuel Delivery Regulation Valve And Injector



**Laboratory Instrument India**