




Email : sales@laboratoryinstrumentindia.com

Product Name :
Common Rail Pump With Injector

Product Code :
LBNY-0005-120000110



Description :

Discharge of the pressurized fuel via a ball valve into the collector line to the rail.
High pressure pump: By turning the eccentric cam the operation of the pump piston can be seen.
Injector: Injection nozzle, solenoid valve for volume control and for control of start of delivery.
Aspiration of the fuel from the interior of the pump via a flat valve.

Technical Specification :

Common Rail Pump With Injector



Laboratory Instrument India