




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Common Rail Pump With Injector

**Product Code :**  
LBNY-0005-120000110



**Description :**

Discharge of the pressurized fuel via a ball valve into the collector line to the rail.  
High pressure pump: By turning the eccentric cam the operation of the pump piston can be seen.  
Injector: Injection nozzle, solenoid valve for volume control and for control of start of delivery.  
Aspiration of the fuel from the interior of the pump via a flat valve.

**Technical Specification :**

Common Rail Pump With Injector



**Laboratory Instrument India**