



Email : sales@laboratoryinstrumentindia.com

Product Name :
Central Door Locking Training Rig

Product Code :
LBNY-0005-120000113



Description :

Designed to understand the operating principles of a central locking system.
Complete with front and rear door lock assemblies and control unit.
Locking, dead-locking and un-locking can be observed during the operation of the Central Locking System Trainer.
Self contained, fully operational Central Locking system trainer.

Technical Specification :

Including door switches, catches, dead locks ignition switch, full set of door locks, boot & fuel cap switches, central locking control unit, interior light & diagnostic socket
Fully operational central door locking trainer manufactured using original components
Based on a GM/OPEL Vectra B
Fitted with manufacturer's diagnostic socket (16 pin)
Pin-out facility on all electrical components using 2mm sockets enabling real-time values to be taken without compromising the wiring harness
Powder coated paint finish
Board mounted in a vertical frame with castors.

Laboratory Instrument India

