



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Central Door Locking Training Rig

**Product Code :**  
LBNY-0005-120000113



**Description :**

Designed to understand the operating principles of a central locking system.  
Complete with front and rear door lock assemblies and control unit.  
Locking, dead-locking and un-locking can be observed during the operation of the Central Locking System Trainer.  
Self contained, fully operational Central Locking system trainer.

**Technical Specification :**

Including door switches, catches, dead locks ignition switch, full set of door locks, boot & fuel cap switches, central locking control unit, interior light & diagnostic socket  
Fully operational central door locking trainer manufactured using original components  
Based on a GM/OPEL Vectra B  
Fitted with manufacturer's diagnostic socket (16 pin)  
Pin-out facility on all electrical components using 2mm sockets enabling real-time values to be taken without compromising the wiring harness  
Powder coated paint finish  
Board mounted in a vertical frame with castors.

---

# Laboratory Instrument India

