



Email : sales@laboratoryinstrumentindia.com

Product Name :
ABS Rig Teves

Product Code :
LBNY-0005-120000117



Description :

The ABS trainer offers an enhanced training resource, hands on training with real components, greater understanding for the student without the constraints of the modern motor vehicle. Fault panel to aid student diagnosis. Diagnostic socket fitted for the connections of fault code readers or blink coding. Pressure gauges show front and rear line pressures and individual pressure at each wheel. Brake pedal kick-back can be observed. Self contained, fully operational ABS trainer, complete with full braking circuit, manufactured using genuine car components. Simulation of wheel lock to activate ABS modulator. Speed simulation 0-160 km/h.

Technical Specification :

Fully operational Anti-Locking Braking System Training Rig, which allows the observation of a working Anti-Locking Brake system via pressure gauges in the brake circuit
Self contained, realistic and be able to simulate, skid or wheel lock with road speeds up to 160 km/h
ABS unit mounted on a trolley stand with castors
Including a Teves Ate system, complete with ABS pump, electronic control unit, wheel sensors, master cylinder and diagnostic socket
Complete with 6 integral fault switches to aid diagnostic testing
Powder coated paint finish
Rotating parts covered by mesh guards.

Dimensions:

Length 1000 mm, Width 1000 mm, Height 1300 mm
Weight: 75 kg.



Laboratory Instrument India