



Email : sales@laboratoryinstrumentindia.com

Product Name :
AC Power Supply Electronic Engineering Lab

Product Code :
LBNY-0002-610018



Description :

AC Power Supply Electronic Engineering Lab

Technical Specification :

AC Power Supply
Compact and light-weight space-saving design
Very quiet due to the adopting naturally-cooled system without a cooling fan
Wide output range
Power factor correction and universal input
Digital display for output voltage and current
Output voltage and current continuously adjustable
Over current protection
Isolative safety
Input voltage: 230 (220) VAC
Output voltage: 0~260 VAC
Protection: over current

AC Power Supply Electronic Engineering Lab, AC Power Supply Electronic Engineering Lab Equipments, AC Power Supply Electronic Engineering Lab Tools, AC Power Supply Electronic Engineering Lab Tool Kits, AC Power Supply Electronic Engineering Lab Manufacturers, AC Power Supply Electronic Engineering Lab Suppliers from India, China, Kenya



Laboratory Instrument India