



Email : sales@laboratoryinstrumentindia.com

Product Name :
Air Conditioning Training Rig R134a

Product Code :
LBNY-0005-120000121



Description :

The Air Conditioning trainer offers an enhanced training resource, hands-on training with real components, greater understanding for the student with out the constraints of the modern motor vehicle.
Self contained fully operational Air Conditioning Trainer.
Complete with compressor, condenser, receiver dryer, expansion valve and evaporator.
Supplied with sight glasses to observe the state of charge of the refrigerant in the high and low pressure sides.
High and low pressure service ports, easily accessed for recovery and re-charging of there frigerant.
Fully portable unit easily transported for off-site training, complete with stand.
Component used represent a modern A/C circuit, environmentally friendly, using R134a refrigerant.

Technical Specification :

Complete with compressor, condenser, receiver dryer, expansion valve and evaporator
Fully operational Air Conditioning Trainer, which allows the observation of a working Air conditioning rig without the constraints of the modern motor vehicle

Refrigerant R134a
Air Conditioning components mounted within a box section frame on a trolley stand with castors
Rotating parts covered by mesh guards
Powder coated paint finish
High and Low pressure service ports
Charge weight 400 grams
Compressor driven by a 240 V 2.2 kW electric motor
Sight glasses fitted, to observe the state of charge of the refrigerant in the High and Low pressure sides.



Laboratory Instrument India