



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Alarm Immobilizer And Central Locking Training Unit

**Product Code :**  
LBNY-0005-120000122



#### **Description :**

Operation of the remote entry, central locking and engine immobilizer system can be observed during the operation of the Alarm and Immobilizer Trainer.

Self contained, fully operational Alarm Immobilizer and Central Locking system trainer.

Designed to understand the operating principles of an alarm and immobilizer system.

Complete with door lock assemblies, alarm and immobilizer modules, key pad entry, high frequency flip or infra-red remote.

Pin-out facility on all sensors and actuators  $\sim$ 2mm socket enabling real-time values to be taken without compromising the wiring harness.

#### **Technical Specification :**

Complete with instrument cluster, interior lights, indicators, siren, atmospheric sensors, door lock assemblies, boot lock, central locking control unit, engine ECU, engine main relay, ignition switch, alarm switch and diagnostic socket

Fully operational Alarm and Immobilizer System Trainer manufactured using original components from RF shield for simulating key transponder error

Powder coated paint finish

Actual door operation using original door latches

Board mounted in a vertical frame with castors.

---

# Laboratory Instrument India

