



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Front Axle Model

**Product Code :**  
LBNY-0005-120000123



**Description :**

The training stand has angle indicator which can be positioned at various locations. Representation of restoring forces. All functions can be changed without tools via eccentric disks, knurled nuts and quick release closures. Front-Axle Model for adjustments with rack and pinion steering.

The following settings are possible:

- Change of the steering trapezoid
- Positive camber
- Negative camber
- Toe-angle
- Toe in
- Toe out
- Follow-up
- Spread

**Technical Specification :**

Dimensions: 167 x 70 x 156 cm, 130 cm  
Weight: 40 kg.

---

# Laboratory Instrument India

