



Email : sales@laboratoryinstrumentindia.com

Product Name :
Airbag Training Model

Product Code :
LBNY-0005-120000124



Description :

This allows working with original manufacturer documents and circuit diagrams. These components are exposed for easy access. Standard diagnostic testers can be connected via the OBD interface. The training model consists of all the original components of a Golf V airbag system. The connections of all functionally relevant sensors and actuators and the control units and components are original connections as in the vehicle. The training equipment contains non-ignitable airbags and belt tensioners. Thus, a safe handling and troubleshooting is possible.

Characteristics:

Original components: Airbag control unit, Driver airbag, Coil spring, Steering column electronics, Passenger airbag, Side airbag, Crash sensor for front and rear side airbag, Instrument cluster with indicator lamps, Belt pretensioner for driver and passenger, Belt buckle for driver and passenger with belt buckle retainer, Passenger seat occupancy recognition (for belt check), Switching off the passenger airbag via key-switch, OBD Diagnosis connection.

Aluminum frame profile on castors with extra-large tabletop, control panels, OBD connector

Fully functional airbag training model, original components VW Golf 5, specially prepared for training purposes.

Technical Specification :

Power supply: 12 V vehicle battery

Dimensions: l x w x h: 1.600 x 930 x 1.660 mm

Weight: 105kg.



Laboratory Instrument India