



Email : sales@laboratoryinstrumentindia.com

Product Name :
Two Section Climatronic

Product Code :
LBNY-0005-120000128



Description :

The structure of the 2-zone Climatronic is based on the current state of the VW Golf. All parts are in original condition. This training system consists of the table model and the corresponding e-learning companion course. The practice course modules deepen and complement the theory modules in a unique way. The unit has a diagnostic port for connection to diagnostic tools. An error box lets you synchronize errors in the sensors and actuators. In addition, numerous test sockets are integrated. The individual control loops can be shown.

Technical Specification :

Two Section Climatronic



Laboratory Instrument India