



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Wheel Alignment Trainer

**Product Code :**  
LBNY-0005-120000129



**Description :**

Adjustable thrust angle  
Adjustable steering angle  
Adjustable camber angle  
Adjustable Toe-in & Toe-out  
Full size or shorter wheel base  
Each vehicle has fully operational and adjustable suspensionsystem, rear axle, steering system, steering wheel, rims and tyres.

**Technical Specification :**

Length 370 cm, Width 175 cm, Height 135 cm  
Weight: 400 kg.



**Laboratory Instrument India**