



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wheel Alignment Trainer

Product Code :
LBNY-0005-120000129



Description :

Adjustable thrust angle
Adjustable steering angle
Adjustable camber angle
Adjustable Toe-in & Toe-out
Full size or shorter wheel base
Each vehicle has fully operational and adjustable suspensionsystem, rear axle, steering system, steering wheel, rims and tyres.

Technical Specification :

Length 370 cm, Width 175 cm, Height 135 cm
Weight: 400 kg.



Laboratory Instrument India