



Email : sales@laboratoryinstrumentindia.com

Product Name :
DC Power Supply Electronic Engineering Lab

Product Code :
LBNY-0002-610019



Description :

DC Power Supply Electronic Engineering Lab

Technical Specification :

DC Power Supply
Output off function
Dual independent adjustable DC outputs.
Independent, series and parallel variety output
Constant voltage and constant current
Low ripple and noise
Four LED displays with 3 digits, selectable for two sets
Auto tracking output.

DC Power Supply Electronic Engineering Lab, DC Power Supply Electronic Engineering Lab Equipments, DC Power Supply Electronic Engineering Lab Tools, DC Power Supply Electronic Engineering Lab Tool Kits, DC Power Supply Electronic Engineering Lab Manufacturers, DC Power Supply Electronic Engineering Lab Suppliers from India, China, Kenya

Laboratory Instrument India

