



Email : sales@laboratoryinstrumentindia.com

Product Name :
Climatronic R134a

Product Code :
LBNY-0005-120000131



Description :

Including evaporator, compressor, condenser, associated pipe work, climatecontrol and electronic control. Fullyoperational climate control trainer manufactured using original components. The system is powered by a 12 V DC battery forthe electrical circuit and by 240 V AC for the speed controlled motor.

Technical Specification :

Based on aCitroen/Peugeot climate control system
Board mountedin a box section frame with castors
Fitted with manufacturers diagnostic socket (16 pin)
Supplyvoltages: 240 V AC and 12 V DC
Powder coatedpaint finish
Supplied witha comprehensive training manual.

Dimensions:
Length 645 mm, Width 1310 mm, Height 1410 mm
Weight: 135 kg

Laboratory Instrument India

