



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Climatronic R134a

**Product Code :**  
LBNY-0005-120000131



**Description :**

Including evaporator, compressor, condenser, associated pipe work, climatecontrol and electronic control. Fullyoperational climate control trainer manufactured using original components. The system is powered by a 12 V DC battery forthe electrical circuit and by 240 V AC for the speed controlled motor.

**Technical Specification :**

Based on aCitroen/Peugeot climate control system  
Board mountedin a box section frame with castors  
Fitted with manufacturers diagnostic socket (16 pin)  
Supplyvoltages: 240 V AC and 12 V DC  
Powder coatedpaint finish  
Supplied witha comprehensive training manual.

Dimensions:  
Length 645 mm, Width 1310 mm, Height 1410 mm  
Weight: 135 kg

---

# Laboratory Instrument India

