



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electric Windows Training Unit

Product Code :
LBNY-0005-120000133



www.laboratoryinstrumentindia.com

Description :

Designed to understand the operating principles of an electrical window system.
Operation of each electric window motor can be observed during the operation of the Electric Window Trainer.
Self contained, fully operational Electric Window System Trainer.
Complete with front and rear electric window assemblies and switches.

Technical Specification :

Length 1000 mm, Width 750 mm, Height 1950 mm
Weight: 50 kg.



Laboratory Instrument India