



Email : sales@laboratoryinstrumentindia.com

Product Name :
Peugeot and Citroen Common Rail Diesel Injection

Product Code :
LBNY-0005-120000135



Description :

Self contained fully operational Petrol and Diesel Engine Rigs. Engines can be supplied with fault insertion boards, this allows a range of open circuit and high resistance faults to be inserted, designed to aid teaching diagnostics, fault finding and troubleshooting.

Complete with start, charge, fuel, air-intake, cooling and exhaust systems.

Petrol and diesel engine rigs are manufactured to replicate the engine bay of the donor vehicle.

The petrol and diesel engine rigs offer an enhanced training resource, hands-on training with real components, greater understanding for the student without the constraints of the modern motor vehicle.

Diagnostic sockets are fitted to the engine rigs as appropriate. Optional Fault Board available.

All rotating parts are fully guarded.

Technical Specification :

Fully operational diesel engine rig, complete with fuel, cooling, turbo and exhaust systems

Engine Rig mounted on a self-contained stand complete with castors

Rotating parts covered by mesh guards

Based on a Turbo Charged Engine, complete with common rail injection, electronic diesel control (CR)

Fitted with manufacturer's diagnostic socket

Powder coated paint finish.

Laboratory Instrument India

