



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
1.9 TDI Turbo Diesel Injection

**Product Code :**  
LBNY-0005-120000137



#### **Description :**

Self contained fully operational Petrol and Diesel Engine Rigs.

The petrol and diesel engine rigs offers an enhanced training resource, hands-on training with real components, greater understanding for the student without the constraints of the modern motor vehicle.

Complete with start, charge, fuel, air-intake, cooling and exhaust systems.

Petrol and diesel engine rigs are manufactured to replicate the engine bay of the donor vehicle.

Diagnostic sockets are fitted to the engine rigs as appropriate. Optional Fault Board available.

All rotating parts are fully guarded.

Engines can be supplied with fault insertion boards, this allows a range of open circuit and high resistance faults to be inserted, designed to aid teaching diagnostics, fault finding and troubleshooting.

#### **Technical Specification :**

Fully operational diesel engine rig, complete with fuel, cooling, turbo and exhaust systems

Fitted with manufacturer's diagnostic socket (16 pin)

Engine Rig mounted on a self-contained stand complete with castors

Powder coated paint finish

Based On APsa/Citroen 2.0 Hdi Turbo Charged Engine, Complete With Common Rail Injection, Bosch Electronic Diesel Control

Rotating parts covered by mesh guards

Dimensions:

Length 1200mm, Width 1200 mm, Height 1300 mm

Weight: 295kg.



**Laboratory Instrument India**