



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Vag Fsi Direct Injection

**Product Code :**  
LBNY-0005-120000138



#### **Description :**

Self contained fully operational Petrol and Diesel EngineRigs.

Engines can be supplied with fault insertion boards, this allows a range of open circuit and high resistance faults to be inserted, designed to aid teaching diagnostics, fault finding and troubleshooting.

Petrol and diesel engine rigs are manufactured to replicate the engine bay of the donor vehicle.

Complete with start, charge, fuel, air-intake, cooling and exhaust systems.

The petrol and diesel engine rigs offers an enhanced training resource, hands-on training with real components, greater understanding for the student with out the constraints of the modern motor vehicle.

Diagnostic sockets are fitted to the engine rigs as appropriate. Optional Fault Board available.

All rotating parts are fully guarded.

#### **Technical Specification :**

Fitted with manufacturers diagnostic socket

Engine Rig mounted on a self-contained stand complete with castors

Fully operational petrol engine rig, complete with fuel, cooling, exhaust and ignition systems

Based on a VAG FSI direct injection

Powder coated paint finish

Rotating parts covered by mesh guards

Dimensions:

Length 1200 mm, Width 1100 mm, Height 1300 mm

Weight: 250 kg.



**Laboratory Instrument India**