



Email : sales@laboratoryinstrumentindia.com

Product Name :
Motronic M1.5 Fuel Injection Dry

Product Code :
LBNY-0005-120000142



Description :

The fuel injector board offers an enhanced training resource, hands-on training with real components, and better understanding for the student without the constraints of the modern motor vehicle.

Self contained, fully operational Fuel Injection Board.

Available in wet running version or dry running, using LED™ to indicate injection operation.

Diagnostic sockets fitted where applicable to the original system.

Pin-out facility on all electrical components using 2 mm sockets enabling real-time values to be taken without compromising the wiring harness.

Complete with all actuators and sensors. Manufactured using components from the original engine build.

Components arranged for ease of use and identification.

Technical Specification :

Motronic M1.5 Fuel Injection Dry



Laboratory Instrument India