



Email : sales@laboratoryinstrumentindia.com

Product Name :
Ignition System Ford DIS

Product Code :
LBNY-0005-120000144



Description :

Self contained, fully operational Ignition Turret trainer, manufactured using real components. With the use of timing light, timing advance and retard can be observed, effected by variable speed and vacuum control. Open circuit and high resistance faults are included on these trainers to aid diagnostic training. Operation of the Ignition System can be observed through the transparent spark plug cover. A variable spark plug gap adjuster will aid in the demonstration of how voltage is affected by an increased air gap. Testing of these Ignition Trainers can be carried out in the same manner as if still fitted to the original vehicle.

Technical Specification :

Ignition System Ford DIS



Laboratory Instrument India