



Email : sales@laboratoryinstrumentindia.com

Product Name :
Petrol Direct Injection CAVD

Product Code :
LBNY-0005-120000148



www.laboratoryinstrumentindia.com

Description :

The engine rig is equipped with a fully functional original engine on the current technological development. The circuit diagrams and repair manuals are included on the supplied CD. Diagnostic sockets are fitted to the engine rigs as appropriate. All rotating parts are fully guarded. The motor is mounted on a movable aluminum frame. The permanently installed test box with 4mm sockets with direct access to all 154 ports of the engine control unit and the adapter cables offer the possibility to represent 20 realistic faults and malfunctions, to locate and fix it. The engine rigs offers an enhanced training resource, hands-on training with real components, greater understanding for the student without the constraints of the modern motor vehicle. In addition, the motor training stand has simulator components (5), an emergency-off switch as well as a table to place a laptop for instance.

Technical Specification :

Fully functional petrol engine, complete with fuel system, cooling and exhaust system
Engine Rig mounted on a self-contained stand complete with castors
Rotating parts covered
Test box with 2mm sockets with direct access to all 154 ports of the ECU
Error box with 20 faults
Original diagnosis socket
Circuit diagrams and repair manuals on CD
Aluminium frame

Dimensions:

Length 1200mm, Width 1100 mm, Height 1300 mm
Weight: 250kg.



Laboratory Instrument India