



Email : sales@laboratoryinstrumentindia.com

Product Name :
Analog Lab Unit Electronic Engineering Lab

Product Code :
LBNY-0002-610021



Description :

Analog Lab Unit Electronic Engineering Lab

Technical Specification :

Analog Lab Unit
Fixed DC: +5V, +12V, -12V
Current Output: 0.3A
Dual DC Power Supply: 0 ~ +25V and 0 ~ -25V
Adjustable AC Power Supply
Voltage Range: 2V-0V-2V, 12V-0V-12V, 15V-0V-15V
Current Output: 500mA
Function Generator:
Output Waveform: Sine, Square, Triangle and TTL
Output Frequency: up to 100KHz in five steps

Analog Lab Unit Electronic Engineering Lab, Analog Lab Unit Electronic Engineering Lab Equipments, Analog Lab Unit Electronic Engineering Lab Tools, Analog Lab Unit Electronic Engineering Lab Tool Kits, Analog Lab Unit Electronic Engineering Lab Manufacturers, Analog Lab Unit Electronic Engineering Lab Suppliers from India, China, Kenya

Laboratory Instrument India

