



Email : sales@laboratoryinstrumentindia.com

Product Name :
Turbo Diesel

Product Code :
LBNY-0005-120000150



Description :

Engines can be supplied with fault insertion boards, this allows a range of open circuit and high resistance faults to be inserted, designed to aid teaching diagnostics, fault finding and troubleshooting.

Self contained fully operational Petrol and Diesel Engine Rigs. Complete with start, charge, fuel, air-intake, cooling and exhaust systems.

Petrol and diesel engine rigs are manufactured to replicate the engine bay of the donor vehicle.

The petrol and diesel engine rigs offers an enhanced training resource, hands-on training with real components, greater understanding for the student without the constraints of the modern motor vehicle.

Diagnostic sockets are fitted to the engine rigs as appropriate. Optional Fault Board available.

All rotating parts are fully guarded.

Technical Specification :

Engine Rig mounted on a self-contained stand complete with castors

Rotating parts covered by mesh guards

Powdercoated paint finish

Fully operational diesel engine rig, complete with fuel, cooling, turbo and exhaust systems

Typically RVA engine code or similar

Based on a Turbo Diesel engine, complete with mechanical diesel injection system

Dimensions:

Length 1200 mm, Width 1200 mm, Height 1300 mm

Weight: 275 kg.



Laboratory Instrument India