



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cylinder Fuel Injection Engine

**Product Code :**  
LBNY-0005-120000152



**Description :**

Engines can be supplied with fault insertion boards, this allows a range of open circuit and high resistance faults to be inserted, designed to aid teaching diagnostics, fault finding and troubleshooting.

Self contained fully operational motorbike engine rigs.

Complete with start, charge, fuel, air-intake, cooling and exhaust systems. The engine rigs are manufactured to replicate the engine bay of the donor bike.

Pin-out facility on all sensors and actuators ~2mm sockets enabling real-time values to be taken without compromising the wiring harness.

**Technical Specification :**

Engine Rig mounted on a self-contained stand complete with castors

Rotating parts covered by mesh guards

Fully operational motorcycle engine rig, complete with fuel, cooling, exhaust and ignition systems

Powder coated paint finish

Based on a two/four cylinder motorcycle engine, complete with fuel injection system.

Dimensions:

Length 1100 mm, Width 1100 mm, Height 1300 mm

Weight: 200 kg.

---

# Laboratory Instrument India

