




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Running Engine

**Product Code :**  
LBNY-0005-120000153



**Description :**

All motors are used engines with low mileage, incl. the circuit diagram and self-study guide. Engine training stand engine OBD, 4-cylinder built on an aluminum frame profile with 4 large castors and all original components, including valves, ignition, wiring harness, controls, immobilizer, diagnostic plug, fuel gauge, fuel pump, complete exhaust system with turbo charger. Rotating or hot engine parts are protected.

**Characteristics:**

Complete with starter, tank, cooling, exhaust, ignition system, OBD diagnosis socket, instrument cluster, emergency stop, etc.

Aluminum frame profile on castors with control panel, rotating and hot parts protected, width approx 900 mm

Schematics, self-study program

Control unit error free

Possibility for exhaust tests

**Technical Specification :**

Power supply: 12 Volt DC / Car battery

Size L x W x H: 169 x 88 x 152 cm

Weight: ca. 310 kg.

---

# Laboratory Instrument India

