




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Exercise Engine

**Product Code :**  
LBNY-0005-120000154



**Description :**

The motors are suitable for timing belt adjustment, timing adjustment, high pressure pump assembly, cylinder head assembly, adjustment of pump-nozzle elements etc. Exercise motor mounted on a mobile frame with castors and brake. At the base there is a oil pan. Many engine types available.

**Technical Specification :**

Exercise Engine



**Laboratory Instrument India**