



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
4 Cylinder Carburetor Engine

**Product Code :**  
LBNY-0005-120000156



**Description :**

Engines can be supplied with fault insertion boards, this allows a range of open circuit and high resistance faults to be inserted, designed to aid teaching diagnostics, fault finding and troubleshooting. Self-contained fully operational motorbike engine rigs. Complete with start, charge, fuel, air-intake, cooling and exhaust systems. Pin-out facility on all sensors and actuators 2mm sockets enabling real-time values to be taken without compromising the wiring harness. The engine rigs are manufactured to replicate the engine bay of the donor bike.

**Technical Specification :**

Engine Rig mounted on a self-contained stand complete with castors  
Rotating parts covered by mesh guards  
Powder coated paint finish  
Fully operational motorcycle engine rig, complete with fuel, cooling, exhaust and ignition systems  
Based on a four cylinder motorcycle engine, complete with carbureted fuel system.



**Laboratory Instrument India**