



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Direct Fuel Injection EFI Proline

**Product Code :**  
LBNY-0005-120000162



#### **Description :**

With this teaching and training rig all functions and processes, errors and measuring values can be demonstrated and examined without emissions as they are in the original engine. The required power supply, wiring and component configuration are already integrated into the system and preset, ready for operation. This means the system can be put into operation immediately upon setup. All original parts and required components are clearly arranged in the IPO-principle on the color printed front panel.

#### **Features:**

- Arrangement in the IPO principle: input - processing - output
- Fault circuit, lockable by a lid
- Can be diagnosed via OBD interface
- Original parts, Original injection pressures
- Measuring of all relevant parts and components in the central circuit diagram measuring field
- Connection for student circuit diagram measuring desks.

#### **Technical Specification :**

Direct Fuel Injection EFI Proline

---

# Laboratory Instrument India

