




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Direct Petrol Injection TSI Proline

**Product Code :**  
LBNY-0005-120000165



#### **Description :**

With this teaching and training rig all functions and processes, errors and measuring values can be demonstrated and examined without emissions as they are in the original engine. The required power supply, wiring and component configuration are already integrated into the system and preset, ready for operation. This means the system can be put into operation immediately upon setup. All original parts and required components are clearly arranged in the IPO-principle on the color printed front panel.

#### **Features:**

- Measuring of all relevant parts and components in the central circuit diagram measuring field
- Fault circuit, lockable by a lid
- Connection for student circuit diagram measuring desks
- Arrangement in the IPO principle: input - processing - output
- Original parts, Original injection pressures
- Can be diagnosed via OBD interface
- Compact device build, mobile, for lab and workshop practice

#### **Technical Specification :**

Dimensions: W x D x H: 141 x 70 x 193 cm  
Power supply: 400 V, 16 A  
Weight: about 180 kg.



**Laboratory Instrument India**