




Email : sales@laboratoryinstrumentindia.com

Product Name :
Training Package of Oscilloscope Sensors Actuators

Product Code :
LBNY-0005-120000169



Description :

This succeeds in bridging the gap from pure theory to the complexity of modern vehicles. Measuring itself is not missed out. Proper evaluation of signals and oscilloscope images allows for the accurate determination of fault causes. Thanks to the diversity of the practical information and the real measured values for the intact and defective vehicle, the present teaching system makes a decisive contribution to the introduction of theory and practical service. So make your trainees to real diagnostic professionals. With the Automotive Sensor Trainer, Automotive Actuator trainer and Automotive Oscilloscope trainer in one package, the helpless replacing of components belongs to the past. In state-of-the-art vehicle systems, on-board-diagnosis is a standard. Workshop practice shows, however, that any input interference recognized by the control unit are only formulated in very vague detail in the error storage. Systematic measuring and the correct interpretation of the measured values are gradually conveyed in extensive educational documents.

Technical Specification :

Power supply: 12 Volt DC
Dimensions: L x W x H 975 x 350 x 100 mm
Weight: About 8 kg



Laboratory Instrument India

