



Email : sales@laboratoryinstrumentindia.com

Product Name :
Training Package Multimeter

Product Code :
LBNY-0005-120000178



www.laboratoryinstrumentindia.com

Description :

The Automotive Electrics Trainer gives step by step basic electrical laws and circuit diagnosis. The use of automotive lamps in the basic circuits and the experience of their non-linear behavior is the first step toward practice. With this compact, versatile learning system students learn the basics of electrics and electronics, electrical basic laws, basic circuits, electronic components and their metrological function check. Even the correct measuring itself is focused as a learning topic. Proper measuring with a multimeter the Automotive Multimeter Trainer creates transparency. Measure the wiring harness, interpret voltages and identify failure causes become a cinch. The Automotive Electronics Trainer provides an infinite variety of information about the metrological behavior of electronic components, on measuring circuits and vehicle-related evaluations.

Technical Specification :

Dimensions: L x W x H 975 x 350 x 100 mm
Power supply: 12 Volt DC
Weight: about 8 kg.



Laboratory Instrument India

