



Email : sales@laboratoryinstrumentindia.com

Product Name :
PLC Starter Kit Electronic Engineering Lab

Product Code :
LBNY-0002-610024



Description :

PLC Starter Kit Electronic Engineering Lab

Technical Specification :

Starter Kit
(CPU with working Memory 1MB for Program, 5MB for Data)
Function Generator
Complete dual channels function/arbitrary waveform generator
Channel independence, coupling, track working mode
200 MSa/s sampling rate and 14-bit vertical resolution per channel
Output of 6 standard waveforms, built-in 50 kinds of arbitrary waveform.

PLC Starter Kit Electronic Engineering Lab, PLC Starter Kit Electronic Engineering Lab Equipments, PLC Starter Kit Electronic Engineering Lab Tools, PLC Starter Kit Electronic Engineering Lab Tool Kits, PLC Starter Kit Electronic Engineering Lab Manufacturers, PLC Starter Kit Electronic Engineering Lab Suppliers from India, China, Kenya



Laboratory Instrument India

