



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Controls System

**Product Code :**  
LBNY-0005-120000180



#### **Description :**

The combination with other system boards enables versatile variants and configuration levels for the control of components which are connected according to the circuit diagram. This module panel represents the center of your T-Varia system in the field of conventional central electrics.

#### **Features:**

Use standard workshop measuring and testing devices for diagnosis and troubleshooting  
85x connectors in 4mm safety design, including 2 x 6mm safety design, for wiring and for measuring  
Connect all original automotive components according to the wiring diagram  
Quick and easy inserting hidden practical errors in the circuit configuration  
Educational documents with tasks and worksheets and associated solutions  
Additional integrated circuit failure box with 20 errors, lockable with lid.

#### **Technical Specification :**

Dimensions: L x H x D: 956 x 250 x 140 mm

Power supply: 12 Volt DC

Weight: about 6 kg.

---

# Laboratory Instrument India

