



Email : sales@laboratoryinstrumentindia.com

Product Name :
Controls System

Product Code :
LBNY-0005-120000180



Description :

The combination with other system boards enables versatile variants and configuration levels for the control of components which are connected according to the circuit diagram. This module panel represents the center of your T-Varia system in the field of conventional central electrics.

Features:

Use standard workshop measuring and testing devices for diagnosis and troubleshooting
85x connectors in 4mm safety design, including 2 x 6mm safety design, for wiring and for measuring
Connect all original automotive components according to the wiring diagram
Quick and easy inserting hidden practical errors in the circuit configuration
Educational documents with tasks and worksheets and associated solutions
Additional integrated circuit failure box with 20 errors, lockable with lid.

Technical Specification :

Dimensions: L x H x D: 956 x 250 x 140 mm
Power supply: 12 Volt DC
Weight: about 6 kg.

Laboratory Instrument India

