



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Radio And Navigation

**Product Code :**  
LBNY-0005-120000182



#### **Description :**

The selection of components is ideal for basic training. Best used to get to know the basic circuits of a motor vehicle infotainment. The radio and navigation system panel is equipped with OEM components from the car.

#### **Components:**

- 1 x Car Rear View Camera with Night Vision
- 1 x Fully automatic 12V motor antenna
- 1 x GPS receiver with long cable
- 1 x Multimedia Navigation car radio with a 18cm / 7 inch touch screen with DVD / CD-ROM drive, USB port, Micro SD, connectors for steering wheel remote control
- 2 x Coaxial speaker 3-way
- 1 x Front panel switch for parking brake
- 1 x Ignition switch with three switch positions.

Antenna connection and status LED's

Connections:

All components can be wired to a total of 37 connector sockets in 4 mm safety design in different colours and all signals can also be measured accordingly. The connections are labeled.

#### **Technical Specification :**

Radio And Navigation



**Laboratory Instrument India**