



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Starter Test Stand

**Product Code :**  
LBNY-0005-120000190



#### **Description :**

The starter can be slowed with a clutch equipped dial to a complete standstill. The resulting currents can be measured via the built-in ammeter. The teaching device starter test stand the principle of a starter can be explained in an impressive way. All necessary components for functioning in the car are installed on the system. The test stand is operated via a conventional 12 V car battery. From these values, the performance of the starter can be determined by calculation.

#### **Components:**

1 x Starter, original VW Polo 9N.  
1 x Shunt 600 A.  
1 x Start button.  
1 x Flywheel with ring gear VW Polo 9N.  
Clutch release bearing complete VW Polo 9N.  
Ammeter to 600 A, after adjustable.

#### **Technical Specification :**

##### **Dimensions and Specifications:**

This table model is mounted on a slotted plate aluminium profile which is fitted with four rubber feet to ensure a secure footing. All structural parts made of coated steel. The system has dimensions of 600 mm x 300 mm x 400 mm. All rotating parts are protected against contact.



**Laboratory Instrument India**